REBOUND^{TM/MC}
heating pad
coussin chauffant

Model / Modèle 970-801-CN

IMPORTANT NOTES

Please read this entire manual carefully before using this device.

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IMPORTANT SAFETY INSTRUCTIONS

DO NOT DESTROY

DANGER: TO REDUCE RISKS OF BURNS, ELECTRIC SHOCK, FIRE, AND PERSONAL INJURY, THIS PRODUCT MUST BE USED IN ACCORDANCE WITH THE FOLLOWING INSTRUCTIONS:

1. READ ALL INSTRUCTIONS.
2. DO NOT USE WHILE SLEEPING.
3. DO NOT USE ON AN INFANT OR ON AN ANIMAL.
4. THIS PAD IS NOT TO BE USED ON OR BY AN INVALID, A SLEEPING OR UNCONSCIOUS PERSON, A PERSON WITH DIABETES, OR A PERSON WITH POOR BLOOD CIRCULATION.
5. DO NOT USE ON AREAS OF INSENSITIVE SKIN.
6. BURNS CAN OCCUR REGARDLESS OF CONTROL SETTING. CHECK SKIN UNDER PAD FREQUENTLY.
7. DO NOT USE IN AN OXYGEN ATMOSPHERE.
8. DO NOT USE PINS OR OTHER METALLIC MEANS TO FASTEN THIS PAD IN PLACE.
9. DO NOT BUNCH, SIT ON OR AGAINST OR CRUSH PAD-AVOID SHARP FOLDS.
10. NEVER PULL THIS PAD BY THE SUPPLY CORD.
11. DO NOT USE THE CORD AS A HANDLE.
12. THIS UNIT SHOULD NOT BE USED BY CHILDREN WITHOUT ADULT SUPERVISION.
13. PLACE PAD ON TOP OF AND NOT UNDER THE PART OF THE BODY NEEDING HEAT.
14. USE THIS PAD ONLY ON A 110-120 VOLT AC CIRCUIT. UNPLUG WHEN NOT IN USE. DO NOT USE WITH GENERATORS, POWER CONVERTERS, OR INVERTERS.
15. DO NOT TAMPER WITH THIS PAD IN ANY WAY.
16. THERE ARE NO USER-SERVICEABLE PARTS. IF FOR ANY REASON THIS PAD DOES NOT FUNCTION SATISFACTORILY, SEE WARRANTY FOR CONSUMER SERVICE INFORMATION.
17. DO NOT USE THIS PAD WITH LINIMENT, SALVE OR OINTMENT PREPARATIONS, ESPECIALLY ONES THAT CONTAIN HEAT-PRODUCING INGREDIENTS. SKIN BURNS COULD RESULT.
18. DO NOT BEND OR PINCH CORD AND DO NOT USE IF CORD IS DAMAGED.
19. NEVER LEAVE THE APPLIANCE UNATTENDED, ESPECIALLY IF CHILDREN ARE PRESENT.
20. LOOP CORD LOOSELY WHEN STORING. TIGHT WRAPPING MAY DAMAGE CORD AND INTERNAL PARTS.
21. WHEN NOT IN USE, STORE THE HEATING PAD IN THE ORIGINAL PACKAGING, IN A DARK, COOL, DRY PLACE.
22. SAVE THESE INSTRUCTIONS.

THIS APPLIANCE HAS A POLARIZED PLUG (ONE BLADE IS WIDER THAN THE OTHER). TO REDUCE THE RISK OF ELECTRIC SHOCK, THIS PLUG IS INTENDED TO FIT IN A POLARIZED OUTLET ONLY ONE WAY. IF THE PLUG DOES NOT FIT FULLY IN THE OUTLET, REVERSE THE PLUG. IF IT STILL DOES NOT FIT, CONTACT A QUALIFIED ELECTRICIAN TO INSTALL THE PROPER OUTLET.

DO NOT ATTEMPT TO MODIFY THIS PLUG OR DEFEAT THIS SAFETY FEATURE IN ANY WAY.


Figure 1


Figure 2


Figure 3

HOW TO USE YOUR XPRESSHEAT™ CONTROLLER

Your Sunbeam Health Rebound™ Heating Pad with XPressHeat™ features a unique fast-heat technology in which you will feel heat in 30 seconds - for faster pain relief!

Make sure Control Connector is securely attached to the Pad Connector (Figure 2). Push both parts together for a tight fit.

NOTE: Control must be connected to heating pad before plugging unit into outlet.

1. After control is connected to pad, plug into 120V AC electrical outlet.
2. To turn the pad on, press the Power Button. Then select from 6 heat settings by scrolling to the desired setting with HEAT▲ or HEAT▼ buttons. The light ring around the controls, called the Heat Meter, will light up one side at the time, indicating progress as the pad heats up. All four sides of the Heat Meter will light up within 30 seconds of powering up, indicating that the heating pad is warm and ready for use.
3. To turn the pad OFF, press the Power Button one more time. The heating pad and lights of the Heat Meter will turn off.
4. DO NOT SIT ON OR AGAINST OR CRUSH PAD – AVOID FOLDING DURING USE.
5. DO NOT USE WHILE SLEEPING.
6. SEE ADDITIONAL WARNINGS.

AUTO-OFF FEATURE

Your heating pad controller includes a selectable Auto-Off feature. By pressing the “Auto Off” button, you can set the timer to “Stay On” continuously or to shut down after 2 hours of continuous use. **NOTE:** After a power outage, the Auto-Off feature automatically resets the timer.

HOW TO USE YOUR REBOUND™ HEATING PAD WITH XPRESSHEAT™

The versatility of your Rebound™ Heating Pad with XPressHeat™ technology allows you to treat most aches and pains with penetrating, soothing heat treatment. The unique weighted layer better conforms to affected areas, providing optimum heat treatment benefits.

1. Place Rebound™ Heating Pad on affected area. Use snaps as necessary to secure into place.
2. Select desired heat setting and auto OFF time on controller.

WASHING INSTRUCTIONS

DO NOT DRY CLEAN! DO NOT BLEACH! DO NOT USE WRINGER FOR DETACHABLE PAD!

1. Disconnect pad from outlet.
2. Detach Pad Connector from the Control Connector.
3. Wash by hand or as follows in an automatic washer using slow agitation and normal spin speed: *Fill washer with warm water. Add minimum amount of all-purpose detergent, agitate to dissolve, measure heating pad before washing. Add heating pad and soak five minutes. Agitate two minutes. Spin. Fill for cool rinse, agitate one minute then spin. Dry in automatic dryer on high heat, five minutes only, or line dry. Remove damp heating pad from dryer, block to original size, and drape to finish drying. Do not iron your heating pad.*
4. Do not plug in pad and pad connector until they are completely dry.

TROUBLESHOOTING	
PROBLEM	SOLUTIONS
<ul style="list-style-type: none">• No Heat• No Power	<ul style="list-style-type: none">• Verify no physical damage to the cord, controller or pad• Verify tight connection between pad and controller.• Verify AC outlet is functional.• Verify cord is properly plugged into AC outlet.• If failure continues, discontinue use. See Warranty Information to contact Consumer Services.
<ul style="list-style-type: none">• Control Shuts Off	<ul style="list-style-type: none">• Auto Off feature is set to turn unit off after selected auto-off time ends.• Re-start unit by pressing On Button.• Verify tight connection between pad and controller.
<ul style="list-style-type: none">• Pad has a Hot Spot	<ul style="list-style-type: none">• Make sure pad is not folded or bunched up in any area.

2 YEAR LIMITED WARRANTY

Sunbeam Products, Inc. doing business as Jarden Consumer Solutions or if in Canada, Sunbeam Corporation (Canada) Limited doing business as Jarden Consumer Solutions (collectively "JCS") warrants that for a period of two years from the date of purchase, this product will be free from defects in material and workmanship. JCS, at its option, will repair or replace this product or any component of the product found to be defective during the warranty period. Replacement will be made with a new or remanufactured product or component. If the product is no longer available, replacement may be made with a similar product of equal or greater value. This is your exclusive warranty. Do NOT attempt to repair or adjust any electrical or mechanical functions on this product. Doing so will void this warranty.

This warranty is valid for the original retail purchaser from the date of initial retail purchase and is not transferable. Keep the original sales receipt. Proof of purchase is required to obtain warranty performance. JCS dealers, service centers, or retail stores selling JCS products do not have the right to alter, modify or any way change the terms and conditions of this warranty.

This warranty does not cover normal wear of parts or damage resulting from any of the following: negligent use or misuse of the product, use on improper voltage or current, use contrary to the operating instructions, disassembly, repair or alteration by anyone other than JCS or an authorized JCS service center. Further, the warranty does not cover: Acts of God, such as fire, flood, hurricanes and tornadoes.

What are the limits on JCS's Liability?

JCS shall not be liable for any incidental or consequential damages caused by the breach of any express, implied or statutory warranty or condition.

Except to the extent prohibited by applicable law, any implied warranty or condition of merchantability or fitness for a particular purpose is limited in duration to the duration of the above warranty.

JCS disclaims all other warranties, conditions or representations, express, implied, statutory or otherwise.

JCS shall not be liable for any damages of any kind resulting from the purchase, use or misuse of, or inability to use the product including incidental, special, consequential or similar damages or loss of profits, or for any breach of contract, fundamental or otherwise, or for any claim brought against purchaser by any other party.

Some provinces, states or jurisdictions do not allow the exclusion or limitation of incidental or consequential damages or limitations on how long an implied warranty lasts, so the above limitations or exclusion may not apply to you.

This warranty gives you specific legal rights, and you may also have other rights that vary from province to province, state to state or jurisdiction to jurisdiction.

How to Obtain Warranty Service

In the U.S.A.

If you have any question regarding this warranty or would like to obtain warranty service, please call 1-888-264-9669 and a convenient service center address will be provided to you.

In Canada

If you have any question regarding this warranty or would like to obtain warranty service, please call 1-800-667-8623 and a convenient service center address will be provided to you.

In the U.S.A., this warranty is offered by Sunbeam Products, Inc. doing business as Jarden Consumer Solutions located in Boca Raton, Florida 33431. In Canada, this warranty is offered by Sunbeam Corporation (Canada) Limited doing business as Jarden Consumer Solutions, located at 20 B Hereford Street, Brampton, Ontario L6Y 0M1. If you have any other problem or claim in connection with this product, please write our Consumer Service Department.

PLEASE DO NOT RETURN THIS PRODUCT TO ANY OF THESE ADDRESSES OR TO THE PLACE OF PURCHASE.

©2015 Sunbeam Products, Inc., doing business as Jarden Consumer Solutions.
All rights reserved.

In the U.S., distributed by Sunbeam Products, Inc., doing business as
Jarden Consumer Solutions, Boca Raton, FL 33431.

In Canada, imported and distributed by Sunbeam Corporation (Canada), Limited,
doing business as Jarden Consumer Solutions, Brampton, ON L6Y 0M1.

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www.sunbeamhealth.ca

For product questions:
Sunbeam Consumer Service
USA: 1.888.264.9669 Ext. 66598
Canada: 1.800.667.8623

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REBOUND^{TM/MC}
heating pad
coussin chauffant

Model / Modèle 970-801-CN

IMPORTANT NOTES

Please read this entire manual carefully before using this device.
Product may vary from picture.

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©2015 Sunbeam Products, Inc., faisant affaire sous le nom de Jarden Consumer Solutions.
Aux États-Unis, distribué par Sunbeam Products, Inc., faisant affaire sous le nom de
Au Canada, importé et distribué par Sunbeam Corporation (Canada), Limited,